

optimis[™] + vitra[™] FUSION multispot laser

Laser 3 in 1
SLT/YAG/Multispot Laser



Glaucoma, Secondary Cataract and Retinal Photocoagulation

optimis[™] + vitra[™] multispot laser

TECHNICAL SPECIFICATIONS





	OPTIMIS FUSION		VITRA MULTISPOT [†]
	SLT	Nd:YAG	Multispot Photocoagulator
Laser source:	Q-switched, frequency doubled Nd:YAG	Q-switched Nd:YAG	Frequency doubled Nd:YAG
Wavelength:	532 nm	1064 nm	532 nm
Energy levels/Power:	0.3 to 2 mJ	0.3 to 10 mJ (1 pulse)	1500 mW (at tissue)
Pulse duration/Exposure:	4 ns	4 ns	10 ms to continuous
Spot Size:	400 microns	10 microns (86.4% energy in diameter) 8 microns FWHM	Monospot: 50 µm to 500 µm Multispot: 100 µm to 500 µm Patterns: square, circle, triple arc, line, single spot
Repetition rate:	2.5 Hz	2 Hz	N/A
Offset: (Anterior and Posterior)	N/A	-150µm, +30µm, +150µm	N/A
Pulse mode:	Single pulse	1, 2 and 3 pulses per shot	Single, repeat, painting, continuous
Aiming beam:	Red 635 nm, adjustable intensity	Red 635 nm, adjustable intensity	635 – 650 nm
Slit lamp:	CSO 980 - LED	CSO 980 - LED	N/A
Magnification:	3 positions: 10x, 16x, 25x	3 positions: 10x, 16x, 25x	N/A
Cooling:	Air cooled	Air cooled	By Peltier effect
Size (laser only):	64.7 (H) x 54 (W) x 48.5 (D) cm 25.5" (H) x 21.3" (W) x 19" (D)	64.7 (H) x 54 (W) x 48.5 (D) cm 25.5" (H) x 21.3" (W) x 19" (D)	14.5 (H) x 18.5 (W) x 30 (D) cm 5.7" (H) x 7.3" (W) x 11.8" (D)
Weight (laser only):	24.15 kg 53.1 lbs	24.15 kg 53.1 lbs	4.7 kg 10.4 lbs
Power requirements:	100-240 VAC, 50/60 Hz	100-240 VAC, 50/60 Hz	100-240 VAC, 250 VA, 50/60 Hz
Delivery Systems:	N/A	N/A	Laser indirect ophthalmoscopes: Heine Omega 500 (integrated laser adapter) Keeler Vantage Plus (external laser adapter)

Specifications are subject to change without notice.

©2015.Quantel Medical, Optimis Fusion and Vitra Multispot are trademarks of Quantel Medical. All rights reserved.

www.quantel-medical.com

Designed and Manufacturer by **Quantel Medical, France**



Headquarters
Quantel Medical
11, rue du Bois Joli - CS40015
63808 Couron d'Auvergne – FRANCE
Tel: +33 (0)4 73 745 745
Fax: +33 (0)4 73 745 700
E-mail: contact@quantel-medical.fr

North America
Quantel USA
601 Haggerty Lane
Bozeman, MT 59715 – USA
Tel: +1 877 782 6835
Fax: +1 406 522 2005
E-mail: info@quantelmedical.com

Representative Offices
Thailand, Chiang Mai
Brazil, Rio De Janeiro

Quantel medical